

PHILIPS

AVENT

Manual Breast Pump

Breastfeeding starter set



SCF430/16



Comfortable pumping, wherever you go

Natural Motion Technology, for quick milk flow

Enjoy the gentle comfort of the Philips Avent portable manual breast pump. With Natural Motion Technology, inspired by baby's own suckling action for a quick letdown. You can easily adjust rhythm and vacuum. Fits practically all nipples.

Simple, gentle and effective

- Helps you release milk fast
- One-size soft adapting silicone cushion

Portable and convenient

- Ideal for moms who pump on-the-go

Comfortable expression

- Easily triggers milk release before or between feeding
- Choose your rhythm, follow your milk flow
- Express without leaning forward
- Express, store and feed

Easy to clean

- Easy to clean and set up

Highlights

Natural Motion Technology



The innovative silicone flange - inspired by baby's unique suckling action - stimulates the nipple to help start your milk flow, fast. For comfortable and effective expression.

Adapts to your nipples



One size fits all. Because we all come in different shapes and sizes, the silicone cushion gently flexes and adapts to fit your unique anatomy. Fits 99.98% of all nipple sizes* (up to 30mm).

Portable and easy to use



The breast pump is small and lightweight, which means it is easy to store and transport, making pumping on-the-go simple and discreet.

Fast milk release



A unique, portable manual breast pump – ideal for mothers who need to gently release their milk while away from home, or before breastfeeding.

Suction and flow control



Feel free to adjust the pumping rhythm and suction to suit your needs, comfort and milk flow.

Sit comfortably



Relax and sit comfortably while you express - thanks to our design that lets you sit upright, rather than having to lean forward.

Easy to clean and set up



Because there are only a few parts, our new pump is easy to set up and clean.

Express, store and feed



Express directly into the storage cup, which prevents leakage or spillage. When baby's hungry, simply attach the soft nipple to the cup for feeding. Save your milk in the fridge or freezer and choose your preferred storage solution.

Specifications

Functions

No leaning forward: Sit in a comfortable position

Soft & adaptive cushion: Gentle stimulation

Portable and lightweight: Easy on the go

Material

Bottle: Polypropylene, BPA free*

Teat: Silicone, BPA free*

Breast pump: BPA free*

What is included

Expression kit with handle: 1 pcs

4oz Bottle with 0m+ nipple: 1 pcs

9oz Bottle with 1m+ teat: 1 pcs

Disposable breast pads: 6 pcs

Sealing disc: 1 pcs

Storage cups and 1x adaptor: 3 pcs

Cushion cover: 1 pcs

Insulation pouch: 1 pcs



* Based on: (1)Mangel et al. Breastfeeding difficulties, breastfeeding duration, maternal body mass index, and breast anatomy: are they related?. Breastfeeding Medicine, 26 April 2019, ((109 participants, Israel); (2)Ziemer et al. Skin changes and pain in the nipple during the 1st week of lactation.

* Journal of Obstetric, Gynecologic & Neonatal Nursing, May 1993, (20 participants (Caucasian), USA); (3)Ramsay et al. Anatomy of the lactating human breast redefined with ultrasound imaging, 2005, (28 participants, Australia).

* BPA Free breast pump: Only associated with the bottle, and other parts that come into contact with breast milk. Following EU regulation, 10/2011